



DNA Appendix G – Instrument/Software Troubleshooting Guide

If you cannot contact your local Power User, please reach out to any Power User for assistance. See the shared google document “Casework Power Users” on the BIO Team Drive.

Contents

Automate Express	2
How to PAUSE a run	2
How to return a tip/go back to starting position	2
Pipette Tip Not Ejected/Jammed	2
Wet Pipette Filter/Replace Pipette Tip	2
No/Partial Initial Sample Lysate Uptake	2
Clogged Mandrels	3
No/Low Elution Volume	4
Hamilton STARlet	4
Vector Error	4
Input file issue during amplification	7
Dropped tip or tip observed on deck	7
Barcode Error	8
Pipetting issues	8
See Power User	10
GeneMapper ID-X	11
Archiving Data	11
Restarting all of the Oracle Services manually	14
File Pathway error	16
STRmix	18
License Error	18



DNA Appendix G – Instrument/Software Troubleshooting Guide

Automate Express

How to PAUSE a run

- 1) Press STOP to pause the run. Open the door, address the issue. Press START to resume the run
 - a) **IMPORTANT!** If you open the instrument door while the instrument is running, the run stops, and it cannot be restarted. If you need to open the instrument door during a run (even in the event the run is stopped due to an error), first press STOP to pause the run, then open the door.

How to return a tip/go back to starting position

- 1) Press STOP.
- 2) Press ESC until you get to the main menu.
- 3) Press 1 for manual.
- 4) Press 2 to return tip. This brings everything back to their starting positions.

Pipette Tip Not Ejected/Jammed

- 1) PAUSE the run
 - a) If run is stopped due to an error, still PAUSE the instrument
 - i) Try to unjam the tip and replace with a new tip
 - b) Press START to resume the run
- 2) Possible causes
 - a) D-ring too much or too little grease
 - b) D-rings need to be replaced

Wet Pipette Filter/Replace Pipette Tip

- 1) Note Instrument and lane of wet filter, notify power user
- 2) Possible causes
 - a) Pipette tip hitting bottom
 - b) Clogged mandrel (use lighted mirror to check)
 - c) Sample volume too low (lower than recommended)


No/Partial Initial Sample Lysate Uptake

- 1) After initial sample uptake step, before any other steps in the process
 - a) PAUSE automate, manual pipette sample into cartridge **well #2**, replace tip for that well, press START to resume the run



DNA Appendix G – Instrument/Software Troubleshooting Guide

b) Cartridge order:



- 2) After run is complete
 - a) No sample uptake
 - i) Notify Power User and TL
 - ii) The analyst is to create a **QIR**
 - (1) Ensure that the workflow includes which AE instrument, tip nozzle number, and sample type (context)
 - iii) Re-extract sample on a different instrument (add appropriate blanks), with TL approval
 - b) Partial sample uptake
 - i) Notify Power user and TL
 - ii) The analyst is to create a **QIR**
 - (1) Ensure that the workflow includes which AE instrument, tip nozzle number, and sample type (context)
 - iii) Depending on the sample type
 - (1) Go forward with quant using the elution recovered, with TL approval
 - (a) Determine if quant value is ok to use for sample type, with TL approval
 - (2) Retain remaining sample volume
 - c) Possible causes
 - i) Pipette tip hitting bottom of tube
 - ii) Wet tip filter during run
 - iii) Clogged mandrel
 - iv) Sample volume too low (lower than recommended)

Clogged Mandrels

- 1) Notify power user
 - a) Use hand held mirror with light to look into the opening of the mandrel
 - i) If brown or white crusty material observed in mandrel opening, note instrument and lane
 - b) Notify DNA analysts to avoid lane for future extraction until the mandrel can be replaced.



DNA Appendix G – Instrument/Software Troubleshooting Guide

No/Low Elution Volume

- 1) No Elution volume
 - a) Check sample tube for remaining sample
 - i) If full sample volume remains, follow “No sample uptake” steps above, with TL approval
 - ii) If partial sample volume remains
 - (1) Depending on the remaining sample volume re-extract, with TL approval
 - (a) If a re-extraction is performed, the analyst is to create a **QIR**
 - (b) Ensure that the workflow includes which AE instrument, tip nozzle number, and sample type (context)
- 2) Low elution volume
 - a) Check sample tube for remaining sample
 - i) Depending on the remaining sample volume re-extract, with TL approval
 - (1) If a re-extraction is performed, the analyst is to create a **QIR**
 - (2) Ensure that the workflow includes which AE instrument, tip nozzle number, and sample type (context)
 - ii) If no sample volume remains, the analyst is to create a **QIR**
 - (1) Ensure that the workflow includes which AE instrument, tip nozzle number, and sample type (context)
- 3) Possible causes
 - a) Initial sample volume too low
 - b) Wet tip filter during run
 - c) Clogged mandrel

Hamilton STARlet

Vector Error

- 1) Before script starts
 - a) Check to make sure there is only one Hamilton window open
 - i) Example error:






DNA Appendix G – Instrument/Software Troubleshooting Guide

b) Check for security updates

i) Example error:



c) Let power user know so they can investigate if the computer was removed from the “do not update” and/or “secops” lists

d) Control panel → programs → view installed updates

i) Look for installed security updates for Microsoft Windows (may start with KB....)

ii) Right click to uninstall one at a time

iii) “Yes” at prompt to uninstall

iv) “Yes” at prompt to restart computer

e) Test connection

i) If vector error received, remove next update

2) During run


a) Immediately after importing sample lists and instructions

i) File name is too long


ii) Example error:



DNA Appendix G – Instrument/Software Troubleshooting Guide




- b) Instrument runs into something unexpected (different tube for example)
- Double check layout of hardware and tubes
 - Did you switch TRIO dilution buffer and standard for example?
 - Example error:



- iv) The Hamilton expected the dilution tube depth, standard tube with shallower hardware definition in its place. Z-axis could not go to prescribed depth as it "smashed." Follow up error after "ok:"




DNA Appendix G – Instrument/Software Troubleshooting Guide



Input file issue during amplification

- 1) May cause Hamilton to only pipette PC and NAB into amp plate
- 2) Ensure that 7500 data export has correct columns chosen and that they are in the correct order, it should look like this



- a) If there are additional columns present or if the column order is different from the above uncheck the boxes on the left side of the screen under “Select Results Content” and re-check them in the correct order prior to exporting quant data

Dropped tip or tip observed on deck

- 1) This may be more common with humidity changes in the laboratory environment
- 2) Procedure to check the deck:




DNA Appendix G – Instrument/Software Troubleshooting Guide

- a) Examine plate
- b) Examine carrier trays & plastics for extraneous liquid
- c) Do a good decontamination of the deck, carrier trays, tube holders etc
- d) Document with a **QIR** to include troubleshooting steps undertaken

Barcode Error

- 1) May occur when loading tip carrier if laser does not properly read the barcodes on the tip racks, or if racks have been placed in carrier backwards (with barcodes facing back of carrier)



- a) Click continue, does not impact run
- b) If happens repeatedly, may indicate that laser needs to be adjusted by a Hamilton service engineer
 - i) Contact the local power user to coordinate

Pipetting issues

- 1) Error messages in dilution table when sample(s) should be processed, but were skipped
 - a) Examples:
 - i) "Sample Concentration higher than 500 ng/ul" and a "-1" value for the total amp target
 - ii) The associated sample file name does not have "_DNU" suffix so your batch notes think it was amped, but the 3500 file generated by the Hamilton skips over it.
 - b) Steps to take:

DNA Appendix G – Instrument/software troubleshooting reference

Document # 60284

Revision # 3

Issue Date: 3/7/2022 11:01:42 AM

Issued by: Lisa Yoshida

Page 8 of 20

All Printed Copies are Uncontrolled



DNA Appendix G – Instrument/Software Troubleshooting Guide

- i) Check input files – are any samples duplicated?
 - ii) Check input files – were edits made to sample names between quant and amp set up?
Check for consistency as 7500 file and amp set up file must match exactly
 - iii) Can ask power user to look with you/second set of eyes
 - iv) Can escalate to TL team if struggling to identify cause
 - v) Otherwise, **QIR** serves as notification to TL team
- 2) Well has incorrect volume visibly evident upon completion of set up:
- a) Examples:
 - i) Master mix present, but no sample/RB (can happen during amp set up)
 - ii) Master mix + some sample/diluent volume (can happen during amp set up)
 - iii) Sample appears present and no master mix
 - b) Steps to take:
 - i) Check input files and dilution table
 - (1) Renamed a sample and not carried through?
 - (2) Duplicate sample or RB names causing issues?
 - (3) Dilution table has a comment indicating sample skipped?
 - ii) If caught before amp:
 - (1) Measure volume present in well and in sample tube after Hamilton set up and provide in **QIR** narrative
 - (2) Add master mix from overage in master mix tube + sample/RB/TE to a new well and update run information
 - (3) Notate in the batch notes what occurred and remediation steps taken
 - c) If caught based on CE data:
 - i) Measure volume present in sample tube after Hamilton set up and provide in **QIR** narrative
 - ii) Reamplify the sample
 - d) Document with a **QIR** including troubleshooting steps performed to ensure casework is not impacted
- 3) Sample is expected to yield a quant value (reference or body fluid indicated) and results in no or low human DNA detected
- a) Re-quant to determine if value is reproducible
 - i) If robust DNA detected, **QIR** for Hamilton issue during original quant set up
 - ii) If low value obtained (consistent with original quant), proceed with amplification, as applicable to the case
 - b) Perform amplification (manual or Hamilton)
 - i) Consider microcon or re-extraction of additional sample
 - ii) Inhibitor could have been co-extracted with DNA
 - iii) Collection issue or inherent to sample (environmentally degraded)

DNA Appendix G – Instrument/software troubleshooting reference

Document # 60284

Revision # 3

Issue Date: 3/7/2022 11:01:42 AM

Issued by: Lisa Yoshida

Page 9 of 20

All Printed Copies are Uncontrolled



DNA Appendix G – Instrument/Software Troubleshooting Guide

- c) Document the circumstances with a **QIR** to include possible root causes and remediation steps taken
- 4) Sample set up (GF/YFP) on Hamilton yields a flat-lined profile that is not consistent with the quantification value
 - a) Check for primer front
 - i) If not present, re-prep the CE plate
 - ii) If present, proceed below
 - b) Check input file(s) and dilution table for correct target/dilution instructions
 - i) Was the Hamilton aspirating <2uL? May be due to current script scheme
 - ii) Dilution table has a comment indicating sample skipped?
 - c) Perform a manual amplification for the desired target of that sample
 - d) Document with a **QIR** including troubleshooting steps performed to ensure casework is not impacted
- 5) Sample set up on Hamilton (GF/YFP) on Hamilton yields a profile of a low quality that is not consistent with the quantification value(s)
 - a) Check input file(s) and dilution table for correct target/dilution instructions
 - i) Was the Hamilton aspirating <2uL? May be due to current script scheme
 - ii) Dilution table has a comment indicating sample skipped?
 - b) Perform a manual amplification for the target of that sample
 - i) Recommend to not use original amplification for replicate analysis in STRmix due to profile quality
 - c) Document with a **QIR** including troubleshooting steps performed to ensure casework is not impacted
- 6) Sample set up (GF/YFP) on Hamilton or manually yields an overblown profile not consistent with a low quantification value(s) obtained from quant set up on Hamilton
 - a) Re-quant (manual or Hamilton) and proceed with subsequent amplification based off new quant value
 - b) Document with a **QIR** reflecting troubleshooting steps taken (including original and new quant values, original and new profiles) to ensure casework is not impacted

See Power User

- 1) Issue in dilution table (unexpected comment, sample being skipped)
 - a) Check input files
 - b) Sample name changes? Duplicate names?
- 2) Final plate does not correspond to expected plate (wells missed/skipped)
 - a) Check input files
 - b) Sample name changes? Duplicate names?
- 3) If cannot figure out, escalate to TL team

DNA Appendix G – Instrument/software troubleshooting reference

Document # 60284

Revision # 3

Issue Date: 3/7/2022 11:01:42 AM

Issued by: Lisa Yoshida

Page 10 of 20

All Printed Copies are Uncontrolled




DNA Appendix G – Instrument/Software Troubleshooting Guide

GeneMapper ID-X

Archiving Data


- 1) First, you will open GMIDX and make sure no projects are open. If you have an open project, save it, and click File>New Project.
- 2) Next, Click Tools>GeneMapper IDX Manager, find your projects (you can sort them alphabetically by holding down "SHIFT" and right clicking the top of the name row).
- 3) For each project that you have, you will export it by clicking the project until it is highlighted and clicking export.



- 4) You will be prompted to find the folder within the IDX folder with the corresponding project name.
 - a) Open that folder and hit save. Do this for all projects.




DNA Appendix G – Instrument/Software Troubleshooting Guide



- 5) You will then delete the project from GMIDX by highlighting only your project and hitting the delete button.
 - a) When you need to re-import a project, you will also use this screen.
 - b) Click on import, find your project in the IDX folder (.ser file) and import it back into IDX after the backup is complete.




DNA Appendix G – Instrument/Software Troubleshooting Guide



- 6) Double check that your .ser files are in the folders.
- 7) Once you have exported all of your .ser files into their corresponding folders, you will go to the IDX folder and move all of your project folders into the ARCHIVE folder at the very top.
 - a) Find your lab and your name and cut and paste your project folders into there.

DNA Appendix G – Instrument/Software Troubleshooting Guide



- b) When Power Users do the backup, they go into the ARCHIVE folder and save everything to an external hard drive.
 - i) Afterwards, they delete all of the projects from that folder.
 - ii) If you need to get data at a later date, you will use the ARCHIVE hard drive to pull that data once a power user has copied it to that folder.

Restarting all of the Oracle Services manually

- 1) Please go to Computer Services
 - a) Control Panel> System and Security (if you are on Win10, this folder may not be present) > Administrative Tools> Services
 - b) Scroll down until you see the Oracle Services
 - i) Example screenshot:




DNA Appendix G – Instrument/Software Troubleshooting Guide

Name	Description	Status	Startup Type	Log On As
Net.Tcp Listen...	Receives ac...		Disabled	Local Servi...
Net.Tcp Port S...	Provides ab...		Disabled	Local Servi...
Netlogon	Maintains a ...	Started	Automatic	Local Syst...
Network Acce...	The Networ...		Manual	Network S...
Network Conn...	Manages o...	Started	Manual	Local Syst...
Network List S...	Identifies th...	Started	Manual	Local Servi...
Network Locat...	Collects an...	Started	Automatic	Network S...
Network Store...	This service ...	Started	Automatic	Local Servi...
Office Source ...	Saves install...		Manual	Local Syst...
Office Softwar...	Enables the ...	Started	Manual	Network S...
Offline Files	The Offline ...	Started	Automatic	Local Syst...
Oracle IFA VSS...		Started	Automatic	Local Syst...
OracleIFAGML...			Manual	Local Syst...
OracleIFAGML...		Started	Automatic	Local Syst...
OracleJobSche...		Started	Automatic	Local Syst...
OracleMTSRec...		Started	Automatic	Local Syst...
OracleServiceL...		Started	Automatic	Local Syst...
Parental Contr...	This service ...		Manual	Local Servi...
Peer Name Re...	Enables ser...		Manual	Local Servi...
Peer Networki...	Enables mul...		Manual	Local Servi...
Peer Networki...	Provides id...		Manual	Local Servi...
Performance C...	Enables re...		Manual	Local Servi...
Performance L...	Performanc...		Manual	Local Servi...
Plug and Play	Enables a c...	Started	Automatic	Local Syst...
PnP-X IP Bus F...	The PnP-X b...	Started	Manual	Local Syst...

- 2) Click on the first Oracle service in the list, and then click Stop
- 3) Then make sure the start-up is set to Automatic, and then click Start:



DNA Appendix G – Instrument/Software Troubleshooting Guide




- Repeat this for each of the Oracle services.
- When done, they should all have "Started" as the Status and "Automatic" for the Start Up type.
- Please note that stopping them and then re-starting is key here
- Now re-launch GMIDX


File Pathway error

- At times projects will open but not appear to be analyzed although all of the edits were previously saved or when reanalyzing it may not recognize the sample.
 - If this occurs GMIDX will need to be reminded of the file pathway.

DNA Appendix G – Instrument/Software Troubleshooting Guide



i) Select Edit -> Define new sample file path



ii) Lead GMIDX through the pathway to the file on the P drive IDX folder where the raw data from the 3500 is saved

iii) This should link the data back to GMIDX




DNA Appendix G – Instrument/Software Troubleshooting Guide


STRmix™

License Error

- 1) When attempting to open STRmix™ and you receive the error: “All licenses are in use or your trial license has expired”
 - a) Close the STRmix™ window.
 - b) Go to the Start menu and select “Floating License Server.”
 - c) Select “Floating License Server Manager.”




- d) If someone else is using STRmix™ on that server, the Floating License Server Manager will show that the license is in use:




1. Click on the license and it will show who is using the license in the Running Clients for Selected License box.



DNA Appendix G – Instrument/Software Troubleshooting Guide



2. If it is someone other than you, that person is either actively using STRmix or did not log off correctly. Check with that person before proceeding.
- e) If no one else is using the server, you will typically see this error:




1. Click "Ok."
2. Select "Connection."
3. Select "Connect to License Server."



DNA Appendix G – Instrument/Software Troubleshooting Guide

- When the window launches, do not change any of the fields. Check the box for “Connect on Startup” if it is not already checked.



- Click “Connect.”
 - The license server will appear in the “Installed License List” section of the Floating License Server Manager.
 - Close the Floating License Server Manager, and relaunch STRmix™.
- f) If you are unable to re-connect the license following the steps above or you continue to get the “Floating License Server Connection Error: Connection Refused,” contact a STRmix™ committee member to have the server restarted.